

Quality | Reliability | Performance

Arthroscopy is an indispensable modality in your practice for diagnostics and therapeutics. Elbows, shoulders, knees/stifles, hips, and more.

When you need the indispensable to be portable, Biovision's **Arthro+ Mobility Suite** has you covered with everything you need in a ready-to-go case:

- Outstanding image performance and capture via the SurgView™ Visualization System, a 21" all-digital imaging platform with integrated touch screen LCD media capture
- Ultra-sharp resolution and brighter-than-ever LED illumination via the HD camera and fiberoptic cable
- Workhorse performance from two arthroscopes – 1.9mm O.D. 60mm W.L. and 2.7mm O.D. 140mm W.L. (with associated cannulas and obturators)
- A full complement of hand tools for grasping/removal, probing, and curettage



Arthro+
Mobility
Veterinary Arthroscopy
ENDOSCOPY SYSTEM





Arthro+ Mobility Veterinary Arthroscopy Endoscopy System Part No. BVTsuite-00009

Item	Description
SurgView™	High definition LED visualization system
Camera	High definition camera with touch membrane-activated controls
MediaCapture USB	Intelligent 21" LCD touch screen monitor w/integrated media capture
Terminal Arm	Terminal/Monitor arm
Fiberoptic Cable	5mm O.D. high performance translucent illumination cable
Coupler	30mm F.L. C-mount coupler
Arthroscope	30 deg., 1.9mm O.D., 60mm W.L.
Cannula	2.78mm O.D., 36.5mm W.L.
Obturator	Conical blunt*
Arthroscope	30 deg., 2.7mm O.D., 140mm W.L.
Cannula	3.97mm O.D., 105mm W.L.
Obturator	Conical blunt*
Forceps	3.5mm, meniscus grasping w/ratchet
Forceps	2.0mm, rigid grasping forceps
Curette	2.0mm spoon
Curette	2.2mm spoon
Probe	90 deg., hook style
Handle	Handle for hook probe
Accessory	Cleaning brush for cannula
Case	Black resin fiber w/casters and telescoping handle
Accessory	Custom foam for SurgView and accessories
Training	Referral to training resources

**Other styles available: pyramidal sharp, spherical blunt*

