



Simply. Better.

Introducing Biovision's Next Generation EndoDiagnostic+Surgical Suite

Experience the evolution of the laparoscopic modality. Biovision's EndoDiagnostic+Surgical Suite (EDSS) proven tool is now simply better, with an all new digital visualization platform featuring our newly developed high-definition camera with more powerful LED illumination and an intelligent LCD touch screen monitor with integrated media capture USB. From basic video otoscopy to advanced laparoscopy, the EDSS suite package scales effortlessly for all of your rigid endoscopy clinical needs.

- HD resolution
- Brighter-than-ever illumination
- Touch screen
- All-digital
- More powerful, useful, and versatile





EndoDiagnostic+Surgical Suite Part No. BVTSUITE-00001HD

Item	Description
SurgView™	High definition LED visualization system
Camera	High definition camera with touch membrane-activated controls
MediaCapture USB	Intelligent LCD touch screen monitor w/integrated media capture
Fiberoptic cable	5mm O.D. high performance translucent illumination cable
SurgAssist™	Insufflator/electrosurgical generator
Laparoscope	11 mm O.D., 6 mm I.D., 270 mm W.L., 2.1 mm channel
Otoscope	0 deg. Oblique, 5.5-20 mm O.D., 90 mm W.L., 2.1 mm channel
Coupler	30 mm F.L. C-mount coupler
EndoBlade (x2)	5.16 mm O.D., 420 mm W.L., bi-polar cautery instrument
Trocar (x2)	11/12 mm, 100 mm W.L.
Trocar	5 mm, 75 mm W.L.
Forceps	5 mm O.D., 430 mm W.L., Babcock w/ratchet
Forceps	5 mm O.D., 430 mm W.L., cup biopsy w/tooth
Handpiece	Monopolar pencil for extracorporeal use
Accessory	Monopolar connecting wire for endoscopic instruments
Hose	Patient insufflation; open end to male Luer lock
Footswitch	Electrosurgical generator activator switch
Regulator	E-tank saddle w/pressure gauge
Hose	CO ² supply; regulator-to-BV quick-connect
Plate	Grounding plate for monopolar cautery
Wires	Wires for grounding plate
Accessory	Cleaning brush for operating laparoscope
Accessory	Cleaning brush for otoscope
DVD	Procedure DVD and client educational materials
Training	Full day of hands-on surgical training

Options

Table	TippyTable™ surgical positioning table
Forceps	5 mm O.D., 330 mm W.L., Babcock w/ratchet
Forceps	10 mm O.D., 330 mm W.L., Babcock w/ratchet
Sterile tray	Clear sterilization tray; 23" x 8" x 3.75"
Sterile drapes	Sterile drapes for camera w/operating laparoscope
Cart	Adjustable cart w/AC power and swivel arm w/21" LCD display

Surgical Training

Hands-on surgical endoscopic training is included with each EndoDiagnostic+Surgical Suite. Some of the most progressive surgical techniques are taught in 8 hours of intensive training at the Center for the Advancement of Rigid Endoscopy (C.A.R.E.) in Denver, Colorado.